

ZIRCONIA ATTACHMENT

TAV advanced and innovative over-denture restoration system offers Zirconia Locker Attachments and Retention Caps allowing high-level functionality, easy-to-use, excellent wear resistance and reliability.



Zirconia attachment connected with titanium base / connector to fit any implant system

	Zirconia Locker NP	Zirconia Locker SP	Zirconia Locker with Ti Base SP	Zirconia Locker Coni. con SP	Zirconia Locker Coni. con WP
Height, mm					
0	7200513	7200290	-	-	-
1	7200514	7200291	-	-	-
2	7200515	7200292	7200608	7201139	7201146
3	7200516	7200293	7200576	7201140	7201147
4	7200517	7200294	7200577	7201141	7201148
5	7200518	7200330	7200578	7201142	7201149
6	7200519	7200331	7200579	7201143	7201150

PROCESSING AND RETENTION COMPONENTS

TAV Retention Caps Set are compatible with Zest Anchors Locator[®]

RETENTION CAPS SET



Low Retention Plastic Cap // Blue // 700gr

4400110



Medium Retention Plastic Cap // Pink // 1400gr

4400111



High Retention Plastic Cap // Transparent // 2200gr

4400112



Black Processing Male Cap

4400116



Spacer Ring

4400117



Denture Cap

7800747



For 20° divergence between implants

EXTENDED RETENTION CAPS SET



Low Retention Plastic Cap // Red // 450gr

4400113



Medium Retention Plastic Cap // Orange // 900gr

4400114



High Retention Plastic Cap // Green // 1800gr

4400115



Black Processing Male Cap

4400116



Spacer Ring

4400117



Denture Cap

7800747



For 40° divergence between implants